



Boğaziçi University

**Introductory
Phys Labs**

1863


REFLECTION AND REFRACTION

PHYL202



REFLECTION AND REFRACTION

In this experiment you will be tracing the light rays reflected or refracted from various optical elements and determine some relevant quantities of these elements.



light source



plane mirror



convex mirror
concave mirror



convex lens



concave lens



prism



1 8 6 3

BOĞAZIÇI ÜNİVERSİTESİ


PART I: REFLECTION

1863

REFLECTION AND REFRACTION

A) PLANE MIRROR


In a plane mirror the incident and reflected angles with respect to the normal are equal



REFLECTION AND REFRACTION


A) PLANE MIRROR

- Please watch the video, and prepare set of white papers, ruler, and protractor to draw your pictures and make your measurements.
- Draw the image on white paper, maximizing your screen lights and turning the lights off in your room will help you.
- With help of protractor, draw the normal line to the surface.
- Measure incident and outgoing angles, θ_i and θ_r , respectively.



REFLECTION AND REFRACTION

Measure Incident and Reflected Angle θ_i and θ_r



A)

Plat

Incident ray angle


$\theta_i = \dots\dots\dots$

Reflected ray angle

$\theta_r = \dots\dots\dots$

REFLECTION AND REFRACTION

Measure Incident and Reflected Angle θ_i and θ_r




- A) **Plan**
- Incident ray angle $\theta_i = \dots\dots\dots$
- Reflected ray angle $\theta_r = \dots\dots\dots$

B) CONCAVE – CONVERGING MIRROR

Focal lengths of concave and convex mirrors are simply half the radius (R) of curvature for the respective surface.

$$f = \frac{R}{2} \quad (\text{R will be determined with help of Chord Method})$$




Arrange size of your video so that 10 cm of the picture matches **exactly** with 10.0 cm of your ruler **as it is already implied in your data video instructions!**

If not possible, you can scale to a smaller value than 10 cm, and rescale your measurements.

B) CONCAVE – CONVERGING MIRROR

Get a white paper and draw the picture while watching your data video

- measure $f_{\text{concave mirror}}$



B) CONCAVE – CONVERGING MIRROR

Fill the empty spaces accordingly.

B) Concave – Converging Mirror:

Focal Length of the mirror f_{EV} =


Radius of the mirror
(From Chord Method) R =

Focal length of the mirror
(From Chord Method) f_{CV} =

% difference in focal lengths = $\frac{|f_{CV} - f_{EV}|}{f_{CV}} \times 100$

C) CONVEX – DIVERGING MIRROR

Focal lengths of concave and convex mirrors are simply half the radius of curvature for the respective surface.




REFLECTION AND REFRACTION

C) CONVEX – DIVERGING MIRROR

Get a white paper and draw the picture

- measure $f_{\text{convex mirror}}$



1863

C) CONVEX – DIVERGING MIRROR

Fill the empty spaces accordingly.

C) Convex – Diverging Mirror:

Focal Length of the mirror f_{EV} =

Radius of the mirror
(From Chord Method) R =

Focal length of the mirror
(From Chord Method) f_{CV} =

Thickness of the mirror x = This is a sanity check, the result should be around ~0.5 cm, otherwise please check your chord method

% difference in focal lengths = $\frac{|f_{CV} - f_{EV}|}{f_{CV}} \times 100$


A large, faint, circular seal of Boğaziçi University is centered in the background. The seal contains the text "BOĞAZIÇI ÜNİVERSİTESİ" around the top and "1863" at the bottom. In the center of the seal is a diamond shape containing a crescent moon and a star.

PART II: REFRACTION

D) CONVEX – CONVERGING LENS

Focal length and the radius of curvature of a lens is related through the following expression:

$$\frac{1}{f} = (n - 1) \frac{2}{R} \quad (\text{R will be determined with help of Chord Method})$$




REFLECTION AND REFRACTION

D) CONVEX – CONVERGING LENS

Get a white paper and draw the picture

- measure $f_{\text{convex lens}}$ *



D) CONVEX – CONVERGING LENS

Fill the empty spaces accordingly.

D) Convex – Converging Lens :

Refraction Index n =
It depends on material. You will read it from the data video.

Focal Length of the lens f_{EV} =

Radius of the convex lens
(From Chord Method) R =


Focal length of the convex lens
(From Chord Method) f_{CV} =

% difference in focal lengths = $\frac{|f_{CV} - f_{EV}|}{f_{CV}} \times 100$

E) CONCAVE - DIVERGING LENS

Focal length and the radius of curvature of a lens is related through the following expression:

$$\frac{1}{f} = (n - 1) \frac{2}{R} \quad (\text{R will be determined with help of Chord Method})$$




REFLECTION AND REFRACTION

E) CONCAVE - DIVERGING LENS

Get a white paper and draw the picture

- measure $f_{\text{concave lens}}$



E) CONCAVE - DIVERGING LENS

Fill the empty spaces accordingly.

E) Concave – Diverging Lens:

Refraction Index n =
 It depends on material. You will read it from the data video.

Focal Length of the lens f_{EV} =



Radius of the concave lens
 (From Chord Method) R =

Focal length of the concave lens
 (From Chord Method) f_{CV} =

% difference in focal lengths = $\frac{|f_{CV} - f_{EV}|}{f_{CV}} \times 100$

CHORD METHOD

- Draw two chords as apart from each other as possible
- Draw perpendicular lines from center of each chord with help of a protractor
- Draw a line from point of intersection of perpendicular lines to the circle



$$f_{CV} = \frac{R}{2}$$

$$\frac{1}{f_{CV}} = (n - 1) \frac{2}{R}$$

$$\% \text{ difference in focal lengths} = \frac{|f_{CV} - f_{EV}|}{f_{CV}} \times 100$$

F) PRISM


We can calculate the index of refraction of a prism by measuring the angle of minimum deviation between incident and outgoing angles, and also prism angle.

$$n = \frac{\sin\left(\frac{D_{min} + A}{2}\right)}{\sin\left(\frac{A}{2}\right)}$$



REFLECTION AND REFRACTION

PRISM




1863

F) PRISM

Get a white paper and draw the picture

Measure minimum deviation angle D_{\min} and prism angle A .



F) PRISM

Fill the empty spaces accordingly.

F) Prism:

Minimum deviation
between incident
and refracted rays

$$D_{\min} = \dots\dots\dots$$

Prism angle

$$A = \dots\dots\dots$$

Index of Refraction

$$n_{EV} = \dots\dots\dots$$

True Value for the
Index of Refraction

$$n_{TV} = \dots\dots\dots$$

% difference for n

$$= \frac{|n_{TV} - n_{EV}|}{n_{TV}} \times 100 \dots$$